

# Streptavidin Fluorescent Conjugates

Caution: For Laboratory Use. A product for research purposes only.

**Product Numbers: NEL750001EA / NEL751001EA**

Streptavidin-HRP Conjugate	NEL750	(2 vials, 0.25mL/vial)
Streptavidin-AP Conjugate	NEL751	(2 vials, 0.25mL/vial)

## STORAGE CONDITIONS

Store at -20°C until first use. Store at 4°C after thawing. **Do not re-freeze this product.**

## DESCRIPTION

Streptavidin-HRP (NEL750) and Streptavidin-AP (NEL751) are highly sensitive detector conjugates. These conjugates bind to biotin labeled proteins in Western blot, ELISA and immunohistochemistry assays. They also bind to nucleic acids labeled with biotinylated nucleotide analogs in Southern and/or Northern blot, colony/plaque screen and in situ hybridization (ISH) assays.

When used in homologous DNA dot hybridizations with a biotin N<sup>6</sup> labeled DNA probe, as little as 100 femtograms of nucleic acid are detectable in 2 hours or CDP-*Star*<sup>®</sup> Nucleic Acid Chemiluminescence Reagent (NEL601) in 5-10 minutes.

## SUGGESTED WORKING CONDITIONS

The following are suggested conjugate concentration ranges for various assays. Individual conditions may require some conjugate titration to obtain optimal signal clarity.

### SA-HRP

Membrane assays (e.g., Southern, Northern, Western, colony/plaque): Dilute conjugate 1:12,500-1:25,000 in 1X Phosphate Buffered Saline (PBS), pH 7.1-7.4 (e.g., Dulbecco's PBS), 0.1% SDS, 0.05% Tween-20, blocking reagent (or equivalent).

Slide assays (i.e., ISH and IHC) including TSA: Dilute conjugate 1:1250-1:2500 in 0.10 M Tris-HCl, pH 7.5, 0.15 M NaCl, blocking reagent (e.g., 2% sheep serum) or equivalent.

### SA-AP

Membrane assays (e.g., Southern, Northern, Western, colony/plaque): Dilute conjugate 1:25,000-1:50,000 in 1X Phosphate Buffered Saline (PBS), pH 7.1-7.4 (e.g., Dulbecco's PBS), 0.1% SDS, 0.05% Tween-20, blocking reagent or equivalent.

Slide assays (i.e., ISH and IHC) without amplification: Dilute conjugate 1:1000-1:2000 in 0.10 M Tris-HCl, pH 7.5, 0.15 M NaCl, blocking reagent (e.g., 2% sheep serum) or equivalent.

Slide assays (i.e., ISH and IHC) with TSA: Dilute conjugate 1:2500-1:5000 in 0.10 M Tris-HCl, pH 7.5, 0.15 M NaCl, blocking reagent (e.g., 2% sheep serum) or equivalent.

## STABILITY

Streptavidin-HRP and Streptavidin-AP are stable for at least 6 months from the date of shipment, when stored at 4°C

## HAZARD WARNING

This product is considered to be non-hazardous. Although the product is classified as non-hazardous, we strongly recommend using prudent laboratory practices: Avoiding unnecessary contact, use of gloves, eye protection, lab coats, etc. when using this or any other laboratory reagent.

CDP-*Star* is a trade mark of Tropix, Inc. CDP-*Star* is protected under one or more of U.S. patents 5,326,882 and 4,931,569.

PerkinElmer, Inc.  
549 Albany Street  
Boston, MA 02118 USA  
P: (800) 762-4000 or (+1) 203-925-4602  
[www.perkinelmer.com/nenradiochemicals](http://www.perkinelmer.com/nenradiochemicals)

**For a complete listing of our global offices, visit [www.perkinelmer.com/ContactUs](http://www.perkinelmer.com/ContactUs)**

Copyright ©2010, PerkinElmer, Inc. All rights reserved.

PerkinElmer® is a registered trademark of PerkinElmer, Inc.

All other trademarks are the property of their respective owners.

