



SARTORIUS ENTRIS BALANCES END OF THE YEAR CAMPAIGN

Entris® II essential line -30% off the price

**BALANCE THAT FITS TO ANY
LABORATORY – ENTRIS® II**

IN SELECTION:

- Weighing capacity 60-8200 g
- Precision 0,1-1000 mg
- Internal calibration
- 12 programs
- RS232 connection
- GLP|GMP printout
- Both verified and regular balances

FEATURES:

- Automatic calibration dependent on time and temperature (isoCAL)
- LED touch-screen
- USB-C connection
- Plug&Play fast and easy connection to PC, no extra software is needed
- BPT surface for high chemical resistance
- Special-modified glass draft shield to minimize weighing errors from static electricity



TAKE A CLOSER LOOK AT ENTRIS® II GENERATION ON SARTORIUS [PAGE](#) AND [BROCHURE](#).

EXAMPLES OF MODELS:

REGULAR MODELS

Article	Name/weighing capacity/precision	Regular price (€)*	Sale price (€)*
<u>BCE64i-1S</u>	Lab balance with int, Cal, 60 g, 0,1 mg	2100.00	1470.00
<u>BCE124i-1S</u>	Lab balance with int, Cal, 120 g; 0,1 mg	2200.00	1540.00
<u>BCE224i-1S</u>	Lab balance with int, Cal, 220 g; 0,1 mg	2330.00	1631.00
<u>BCE223i-1S</u>	Lab balance with int, Cal, 220 g; 1 mg	1400.00	980.00
<u>BCE323i-1S</u>	Lab balance with int, Cal, 320 g; 1 mg	1520.00	1064.00
<u>BCE423i-1S</u>	Lab balance with int, Cal, 420 g; 1 mg	1630.00	1141.00
<u>BCE623i-1S</u>	Lab balance with int, Cal, 620 g; 1 mg	1860.00	1302.00
<u>BCE653i-1S</u>	Lab balance with int, Cal, 650 g; 1 mg	2100.00	1470.00
<u>BCE622i-1S</u>	Lab balance with int, Cal, 620 g; 10 mg	1270.00	889.00
<u>BCE822i-1S</u>	Lab balance with int, Cal, 820 g; 10 mg	1400.00	980.00
<u>BCE1202i-1S</u>	Lab balance with int, Cal, 1200 g; 10 mg	1520.00	1064.00
<u>BCE2202i-1S</u>	Lab balance with int, Cal, 2200 g; 10 mg	1630.00	1141.00
<u>BCE3202i-1S</u>	Lab balance with int, Cal, 3200 g; 10 mg	1750.00	1225.00
<u>BCE4202i-1S</u>	Lab balance with int, Cal, 4200 g; 10 mg	1860.00	1302.00
<u>BCE6202i-1S</u>	Lab balance with int, Cal, 6200 g; 10 mg	2100.00	1470.00
<u>BCE2201i-1S</u>	Lab balance with int, Cal, 2200 g; 100 mg	1400.00	980.00
<u>BCE5201i-1S</u>	Lab balance with int, Cal, 5200 g; 100 mg	1520.00	1064.00
<u>BCE8201i-1S</u>	Lab balance with int, Cal, 8200 g; 100 mg	1630.00	1141.00
<u>BCE6200i-1S</u>	Lab balance with int, Cal, 6200 g; 1 g	1400.00	980.00
<u>BCE8200i-1S</u>	Lab balance with int, Cal, 8200 g; 1 g	1520.00	1064.00

*VAT will be added.
Campaign prices
are valid
until 31.12.2023.



VERIFIED MODELS

Article	Name/weighing capacity/precision	Regular price (€)*	Sale price (€)*
<u>BCE64i-1CEU</u>	Lab balance with int, Cal, 60 g, 0,1 mg	2570.00	1799.00
<u>BCE124i-1CEU</u>	Lab balance with int, Cal, 120 g; 0,1 mg	2690.00	1883.00
<u>BCE224i-1CEU</u>	Lab balance with int, Cal, 220 g; 0,1 mg	2800.00	1960.00
<u>BCE223i-1CEU</u>	Lab balance with int, Cal, 220 g; 1 mg	1860.00	1302.00
<u>BCE323i-1CEU</u>	Lab balance with int, Cal, 320 g; 1 mg	1980.00	1386.00
<u>BCE423i-1CEU</u>	Lab balance with int, Cal, 420 g; 1 mg	2100.00	1470.00
<u>BCE623i-1CEU</u>	Lab balance with int, Cal, 620 g; 1 mg	2330.00	1631.00
<u>BCE653i-1CEU</u>	Lab balance with int, Cal, 650 g; 1 mg	2570.00	1799.00
<u>BCE622i-1CEU</u>	Lab balance with int, Cal, 620 g; 10 mg	1750.00	1225.00
<u>BCE822i-1CEU</u>	Lab balance with int, Cal, 820 g; 10 mg	1860.00	1302.00
<u>BCE1202i-1CEU</u>	Lab balance with int, Cal, 1200 g; 10 mg	1980.00	1386.00
<u>BCE2202i-1CEU</u>	Lab balance with int, Cal, 2200 g; 10 mg	2100.00	1470.00
<u>BCE3202i-1CEU</u>	Lab balance with int, Cal, 3200 g; 10 mg	2220.00	1554.00
<u>BCE4202i-1CEU</u>	Lab balance with int, Cal, 4200 g; 10 mg	2330.00	1631.00
<u>BCE6202i-1CEU</u>	Lab balance with int, Cal, 6200 g; 10 mg	2570.00	1799.00
<u>BCE2201i-1CEU</u>	Lab balance with int, Cal, 2200 g; 100 mg	1860.00	1302.00
<u>BCE5201i-1CEU</u>	Lab balance with int, Cal, 5200 g; 100 mg	1980.00	1386.00
<u>BCE8201i-1CEU</u>	Lab balance with int, Cal, 8200 g; 100 mg	2100.00	1470.00
<u>BCE6200i-1CEU</u>	Lab balance with int, Cal, 6200 g; 1 g	1860.00	1302.00
<u>BCE8200i-1CEU</u>	Lab balance with int, Cal, 8200 g; 1 g	1980.00	1386.00

*VAT will be added.

Campaign prices are valid until 31.12.2023.

Ask more from our
product manager:

Karolis Šalčiūnas

+370 61718204

karolis@novabio.lt



SARTORIUS

Authorized Partner

Sartorius and the Sartorius logo are registered trademarks of Sartorius AG or any of its affiliates.