



# SARTORIUS BALANCES END OF THE YEAR CAMPAIGN

## Entris® II & Cubis® II -30% off the price

BALANCE THAT FITS TO ANY  
LABORATORY – NEW **ENTRIS® II**

### IN SELECTION:

- Weighing capacity 60-8200 g
- Precision 0,1-1000 mg
- Internal calibration
- 12 programs
- RS232 connection
- GLP|GMP printout
- Both verified and regular balances

### WHAT IS NEW?

- Automatic calibration dependent on time and temperature(isoCAL)
- LED touchscreen
- USB-C connection
- Plug&Play fast and easy connection to PC, no extra software is needed
- BPT surface for high chemical resistance
- Special-modified glass draft shield, to minimize weighing errors from static electricity



TAKE A CLOSER LOOK AT NEW ENTRIS® II GENERATION ON SARTORIUS [PAGE](#) AND [BROCHURE](#).

## EXAMPLES OF MODELS:

### REGULAR MODELS

Article	Name/weighing capacity/precision	Regular price (€)*	Sale price (€)*
<b><u>BCE64i-1S</u></b>	Lab balance with int, Cal, 60 g, 0,1 mg	2000.00	1400.00
<b><u>BCE124i-1S</u></b>	Lab balance with int, Cal, 120 g; 0,1 mg	2110.00	1477.00
<b><u>BCE224i-1S</u></b>	Lab balance with int, Cal, 220 g; 0,1 mg	2220.00	1554.00
<b><u>BCE223i-1S</u></b>	Lab balance with int, Cal, 220 g; 1 mg	1330.00	931.00
<b><u>BCE323i-1S</u></b>	Lab balance with int, Cal, 320 g; 1 mg	1450.00	1015.00
<b><u>BCE423i-1S</u></b>	Lab balance with int, Cal, 420 g; 1 mg	1550.00	1085.00
<b><u>BCE623i-1S</u></b>	Lab balance with int, Cal, 620 g; 1 mg	1770.00	1239.00
<b><u>BCE622i-1S</u></b>	Lab balance with int, Cal, 620 g; 10 mg	1210.00	847.00
<b><u>BCE822i-1S</u></b>	Lab balance with int, Cal, 820 g; 10 mg	1330.00	931.00
<b><u>BCE1202i-1S</u></b>	Lab balance with int, Cal, 1200 g; 10 mg	1450.00	1015.00
<b><u>BCE2202i-1S</u></b>	Lab balance with int, Cal, 2200 g; 10 mg	1550.00	1085.00
<b><u>BCE3202i-1S</u></b>	Lab balance with int, Cal, 3200 g; 10 mg	1670.00	1169.00
<b><u>BCE4202i-1S</u></b>	Lab balance with int, Cal, 4200 g; 10 mg	1770.00	1239.00
<b><u>BCE6202i-1S</u></b>	Lab balance with int, Cal, 6200 g; 10 mg	2000.00	1400.00
<b><u>BCE2201i-1S</u></b>	Lab balance with int, Cal, 2200 g; 100 mg	1330.00	931.00
<b><u>BCE5201i-1S</u></b>	Lab balance with int, Cal, 5200 g; 100 mg	1450.00	1015.00
<b><u>BCE8201i-1S</u></b>	Lab balance with int, Cal, 8200 g; 100 mg	1550.00	1085.00
<b><u>BCE6200i-1S</u></b>	Lab balance with int, Cal, 6200 g; 1 g	1330.00	931.00
<b><u>BCE8200i-1S</u></b>	Lab balance with int, Cal, 8200 g; 1 g	1450.00	1015.00

\*VAT will be added.  
Campaign prices  
are valid  
until 31.12.2022.

## VERIFIED MODELS

Article	Name/weighing capacity/precision	Regular price (€)*	Sale price (€)*
<u><a href="#">BCE64i-1CEU</a></u>	Lab balance with int, Cal, 60 g, 0,1 mg	2450.00	1715.00
<u><a href="#">BCE124i-1CEU</a></u>	Lab balance with int, Cal, 120 g; 0,1 mg	2560.00	1792.00
<u><a href="#">BCE224i-1CEU</a></u>	Lab balance with int, Cal, 220 g; 0,1 mg	2670.00	1869.00
<u><a href="#">BCE223i-1CEU</a></u>	Lab balance with int, Cal, 220 g; 1 mg	1770.00	1239.00
<u><a href="#">BCE323i-1CEU</a></u>	Lab balance with int, Cal, 320 g; 1 mg	1890.00	1323.00
<u><a href="#">BCE423i-1CEU</a></u>	Lab balance with int, Cal, 420 g; 1 mg	2000.00	1400.00
<u><a href="#">BCE623i-1CEU</a></u>	Lab balance with int, Cal, 620 g; 1 mg	2220.00	1554.00
<u><a href="#">BCE622i-1CEU</a></u>	Lab balance with int, Cal, 620 g; 10 mg	1670.00	1169.00
<u><a href="#">BCE822i-1CEU</a></u>	Lab balance with int, Cal, 820 g; 10 mg	1770.00	1239.00
<u><a href="#">BCE1202i-1CEU</a></u>	Lab balance with int, Cal, 1200 g; 10 mg	1890.00	1323.00
<u><a href="#">BCE2202i-1CEU</a></u>	Lab balance with int, Cal, 2200 g; 10 mg	2000.00	1400.00
<u><a href="#">BCE3202i-1CEU</a></u>	Lab balance with int, Cal, 3200 g; 10 mg	2110.00	1477.00
<u><a href="#">BCE4202i-1CEU</a></u>	Lab balance with int, Cal, 4200 g; 10 mg	2220.00	1554.00
<u><a href="#">BCE6202i-1CEU</a></u>	Lab balance with int, Cal, 6200 g; 10 mg	2450.00	1715.00
<u><a href="#">BCE2201i-1CEU</a></u>	Lab balance with int, Cal, 2200 g; 100 mg	1770.00	1239.00
<u><a href="#">BCE5201i-1CEU</a></u>	Lab balance with int, Cal, 5200 g; 100 mg	1890.00	1323.00
<u><a href="#">BCE8201i-1CEU</a></u>	Lab balance with int, Cal, 8200 g; 100 mg	2000.00	1400.00
<u><a href="#">BCE6200i-1CEU</a></u>	Lab balance with int, Cal, 6200 g; 1 g	1770.00	1239.00
<u><a href="#">BCE8200i-1CEU</a></u>	Lab balance with int, Cal, 8200 g; 1 g	1890.00	1323.00

\*VAT will be added.  
Campaign prices  
are valid  
until 31.12.2022.

## NEW CUBIS® II - THE NEW GENERATION OF MODULAR PREMIUM BALANCES

### IN SELECTION:

- Choose from a range of 45 different weighing modules that best fit your preferences.
- Choose from 7 different draft shields options, depending on your balance type – each with unique features to support maximal balance performance and user-centric design.

### EXAMPLES OF MODELS:

Article	Regular price (€)*	Sale price (€)*
<b>MCE10202S-2S00-O</b>	<b>7402,00</b>	<b>5181,40</b>
MCE panel   Weighing module 1202S, max load 10200 g, precision 0,01 g   Manual leveling Standard model   Without draft shield   Standard application package (10 applications)		
<b>MCA2.7S-2CEU-M-QP1</b>	<b>37 401,00</b>	<b>26 180,70</b>
MCA panel   Weighing module 2.7S, max load 2,1 g, precision 0,0000001 g   Full automatic leveling EU approval, verified   Automatic motorized draft shield   Pharma application package		

\*VAT will be added. Campaign prices are valid until 31.12.2022.



Ask more from our  
product manager:

**Sten Sild**

+372 5367 7740

[sten@novabio.ee](mailto:sten@novabio.ee)

TAKE A CLOSER LOOK AT NEW CUBIS® II GENERATION ON  
SARTORIUS [PAGE](#) AND CONFIGURE YOUR BALANCE [HERE](#).

**SARTORIUS**

Authorized Partner

Sartorius and the Sartorius logo are registered  
trademarks of Sartorius AG or any of its affiliates.