

Technical Datasheet

		EN ISO 13485 certified
Article - No:	07-369-2015	Revision: 4
Remarks		single use only

Specification

Product	Pipette tips PP, premium surface, 100-1.250µl, 96 pcs/rack, transparent, moulded rings, np pcr ready, autoclavable up to 121°C, CE/IVD
Classification	CE/IVD
Dimensions	length: 98,40 mm (+/- 0,2 mm) inner diameter: 0,89 mm (+/- 0,2 mm)
Colour	transparent
Graduation	moulded rings
Nominal Capacity	100-1.250µl
Material Characteristics	polypropylene (PP), premium surface
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62 sterile SAL 10-3
Purity Grade	free of RNase and Dnase free of pyrogens free of PCR inhibitors free of human DNA free of endotoxins free of ATP
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 768 pcs per carton: 3.840 pcs per pallet: 107.520
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C

Technical Datasheet

Date: 18.02.20

Page 2 of 2

pal. not stackable

Approved by nerbe plus quality department
--