

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>05-061-0000</b>	<b>Revision: 2</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Analyzer cup PS, 1,5ml, Ø13,8x22,6 mm, for autoanalyzer from Technikon, Gemsac, Kone and Astra, transparent, CE/IVD
<b>Classification</b>	CE/IVD
<b>Dimensions</b>	diameter: 13,80 mm (+/- 0,2 mm) height: 22,60 mm (+/- 0,2 mm)
<b>Colour</b>	transparent
<b>Nominal Capacity</b>	1,5ml
<b>Material Characteristics</b>	polystyrene (PS)
<b>Product Characteristics</b>	biologically inert physiologically harmless no heavy metals according to EC/94/62 Food contact: raw material in compliance with regulation EC 1935/2004
<b>Shelf - Life</b>	60 months at controlled storage conditions
<b>Packaging</b>	pcs per bag: 1.000 pcs per carton: 10.000 pcs per pallet: 120.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
<b>Approved by nerbe plus quality department</b>	