

Technical Datasheet

EN ISO 13485 certified	
Article - No:	05-181-0000
Remarks	Revision: 2 single use only

Specification

Product	Analyzer cup PS, 3ml, Ø16,75x38 mm, for Hitachi and Boehringer (705, 706, 712,7150, 7250), transparent
Dimensions	diameter: 16,75 mm (+/- 0,2 mm)
	height: 38,00 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	3ml
Material Characteristics	polystyrene (PS)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 1.000
	pcs per carton: 6.000
	pcs per pallet: 144.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

