

Technical Datasheet

Date: 22.06.20

Page 1 of 1

EN ISO 13485 certified		
Article - No:	02-351-6052	Revision: 5
Remarks		single use only

Specification

Product	Test tube PS with fitted screw cap PE, 12ml, Ø16x100 mm, round bottom, tube transparent, cap red, sterile R, max. RCF 3.000g
Dimensions	outer diameter: 16,00 mm (+/- 0,2 mm) height: 100,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap red
Nominal Capacity	12ml
Material Characteristics	polystyrene (PS)
Temperature range	-10°C – 70°C
Centrifugation	max. RCF 3.000g
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Purity Grade	do not contain any pyrogens
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 200 pcs per carton: 1.000 pcs per pallet: 24.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	