

# Technical Datasheet

|                        |                        |
|------------------------|------------------------|
| EN ISO 13485 certified |                        |
| <b>Article - No:</b>   | <b>11-802-1150</b>     |
| <b>Revision:</b>       | <b>1</b>               |
| <b>Remarks</b>         | <b>single use only</b> |

## Specification

|  |  |
|--|--|
| <b>Product</b>                                   | Multi purpose container PP, 200ml, Ø66x88 mm, with enclosed snap-on lid PE, transparent, graduated |
| <b>Dimensions</b>                                | outer diameter: 66,00 mm (+/- 0,2 mm)  |
|  | inner diameter: 62,00 mm (+/- 0,2 mm)  |
|  | height (with lid): 88,00 mm (+/- 0,2 mm)   |
|  | height (without lid): 87,00 mm (+/- 0,2 mm)  |
| <b>Colour</b>                                    | transparent  |
| <b>Graduation</b>                                | grad.  |
| <b>Nominal Capacity</b>                          | 200ml  |
| <b>Material Characteristics</b>                  | container: polypropylene (PP)  |
|  | cap: polyethylene (PE)   |
| <b>Temperature range</b>                         | -80°C - 75°C   |
| <b>Product Characteristics</b>                   | biol. inert  |
|  | physiol. harmless  |
|  | no heavy metals acc. to EC/94/62   |
|  | Food contact: raw material in compl. with reg. EC 1935/2004  |
| <b>Shelf - Life</b>                              | 25 years at controlled storage conditions  |
| <b>Packaging</b>                                 | pcs per bag: 660   |
|  | pcs per carton: 660  |
|  | pcs per pallet: 10.560   |
| <b>Storage</b>                                   | dry condition: 40%-60% rel. humidity   |
|  | bags / cases orig. closed  |
|  | storage temp. +10°C to max. +30°C  |
|  | pal. not stackable   |
| <b>Approved by nerbe plus quality department</b> |  |

00-000-0001/06.2020

