

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>05-900-1050</b>	<b>Revision: 2</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Cryo vial PP, screw cap PE, 5ml, self-standing, transparent, grad., sterile R, autocl.
<b>Dimensions</b>	height: 88,5mm
<b>Colour</b>	transparent
<b>Graduation</b>	grad.
<b>Nominal Capacity</b>	5ml
<b>Material Characteristics</b>	tube: polypropylene (PP) cap: polyethylene (PE)
<b>Temperature range</b>	-196°C - 121°C
<b>Autoclavability</b>	autocl.
<b>Purity Grade</b>	do not contain any pyrogens
<b>Sterilisation Process</b>	radiation sterilization SAL 10-6
<b>Packaging</b>	pcs per bag: 100 pcs per carton: 500 pcs per pallet: 49.500
<b>Storage</b>	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable
<b>Approved by nerbe plus quality department</b>	

00-000-0001/06.2020

