

Phosphatase Labeled Antibodies

Research Use Only. Not for use in diagnostic procedures.

Phosphatase Labeled Antibody to Mouse IgG (H+L) (HAS) [Goat]

Product Number: NEF824001EA

SIZE: 1.0 ml

PRODUCT DESCRIPTION: Affinity purified antibody, isolated from a pool of serum goats immunized with purified Mouse IgG was labeled with phosphatase using a modified glutaraldehyde procedure.

FORM/STORAGE: Liquid. Store @ 2- 8°C

CONTENTS: Each vial contains 1.0 ml @ 1.0 mg/ml of affinity purified Goat Anti-Mouse IgG .

E/P RATIO: Molar enzyme / antibody protein ration = 2.7 : 1.

SPECIFICITY: Tested by immunoelectrophoresis, gel diffusion & ELISA techniques as applicable. This product reacts specifically with Mouse IgG and may recognize other immunoglobulin types that have light chains in common with IgG. Cross-reactivity to human serum has been minimized using affinity procedures. No antibody is detectable to non-immunoglobulin serum components.

STABILIZERS: Stabilization diluents added as a protein stabilizer. No preservatives added.

STABILITY: This product is stable for one year from date of receipt when stored properly.

PerkinElmer, Inc.
549 Albany Street
Boston, MA 02118 USA
P: (800) 762-4000 or (+1) 203-925-4602
www.perkinelmer.com/nenradiochemicals
For a complete listing of our global offices, visit www.perkinelmer.com/ContactUs
Copyright ©2010, PerkinElmer, Inc. All rights reserved.
PerkinElmer® is a registered trademark of PerkinElmer, Inc.
All other trademarks are the property of their respective owners.

